



Waxa la sameeyo haddii...

Haddii aarmi dhaco:

 Aamusii alaarmiga (haddii ay daruuri tahay)


Baadh sababta dabeena ku-soco ficilka habboon

 Sii-daa alarmiga  Ku-wad quudinta




Baateri madhan

Islamarkiiba ku-xidh bamka goobta wyn ee qaybisada.




Faarqi kiishka/hawada ee habaynta alaarmiga

Beddel weelka oo habee ama ka hubi xunbooyin hawo ah.




Bilow Xusuusiyaha

Ku-soco tallabada xigta ama dammi bamka.




Alaarmiga joojinta koru-socodka

Hubi habaynta bamka /tubadda qiimeynta siinta inay saafi tahay kadib bamka, ee bukaanka garabkiisa.



Alaarmiga furitaanka albaabka

Xidh albaabka bamka.




Alaarmiga joojinta koru-socodka

Hubi in habbaynta alaarmiga saafi yahay inta udhaxaysa weelka iyo bamka.




Fariinta ku-dhawaad gaarista mugga

Sii-wad ama joojinta siinta.




Habayntarakibaadda alaarmiga

Hubi goobta habbaynta ka kore iyo hoose nidaamka bamka oo si sax ah u geli.



Mugga bartilmaameedka ayaa gaaray alaarmiga

Dhammaadka quudinta.



Alaarmiga qaladka bamka

Qor koodhka qaladka teknikada ah (Err xyz) oo laxiriir waaaxda daawaynta-da shaqsi.



**FRESENIUS
KABI**

caring for life

Amika®

ENTERAL NUTRITION PUMP HAGAHA TIXRAACA SAWIRKA



WIXII SU'AALO AMA MACLUUMAAD AH WAC HOMECARE 24 SAAC
LAMBARKA KHADKA TELEFOONKA 0808 100 1990.



**FRESENIUS
KABI**

caring for life

Fresenius Kabi Limited,
Cestrian Court, Eastgate Way,
Manor Park, Runcorn,
Cheshire WA7 1NT
telephone: 01928 533533 fakis: 01928 533520
email: scientific.affairsUK@fresenius-kabi.com
www.fresenius-kabi.com/gb

Taariikhda diyaarinta: Marso 2018 EN01046

Amika® waa caalaamad ganacsi oo diwaangashan ee
Fresenius Kabi AG.
Fresenius Kabi Ltd waa isticmaalaha loo oggol yahay ©
Fresenius Kabi Ltd.

01



Isku xir qalaba celisima eek u jira sanduuqa lagu siiyay kadib Ku xir EasyBag ama weel.

02



Fur albaabka is celiya.

03



Gali qalabka (pinch clamp) laynka hoose ilaat"ay guriso" goobta.

04



Ku dheji tuubo bixinta tooska ah gudaha tubada hagaha h.

05



Albaabka xir, hoos u riix si aad u xirto meel adag.

06



Shid Amika® Bangaree - adigoo riixaya Botonka shidista/daminta.

07



Riixfuraha atomaatiga ah, oo kadib riix badhanka bilaw.

08



Marka hawshu dhammaato, riix mar labaad badhanka daarista.

09



Samee xaddiga socodka aad rabto adiga oo riixaya furayaasha +/- illaa inta mugga loo baahan yahay la gaarayo.

10



Samee mugga bartilmaameedka ah adigoo riixaya furayaasha +/- ilaa mugga la doonayo la gaaro.

11



Riix furo kicina si u daarto bilawdo quudinta.

12



Si aad u joojiso quudinta waqti kasta, riix badhanka joogsiga. Marka quudinta dhammaato ama mugga bartilmaameedka la gaaray, riix badhanka jooji.

13



Haddii aad u baahan tahay inaad dib u bilowdo barashada horumarka ka dib marka mugga bartilmaameedka la gaaro, riix badhanka furaha 2-3 ilbiriqsi.

14



Si aad u demiso mishiinka badhka dabka / korontada ilaa inta dib loo soo celinayo.

