15 mm 138 mm 3 mm

10 L NaCl 0,9% Fresenius

UK - NaCl 0.9 % Fresenius, Solution for irrigation
9 mg/ml Sodium Chloride
Air-vent needle behind booklet.

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Seite 1

15 mm 138 mm 3 mm

Close the roller clamp of transfer set if applicable. Insert the spike centered and straight until red rubber stopper meets the shoulder of the spike.

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Open the package of the venting needle. The venting needle can be found in the compartment behind the label. Remove the cap. Do not touch the needle. Insert the needle into the venting point on the hanging/lifted Kabi Container (indicated by arrows) centered and straight all the way.

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Leave the venting needle all the way in the Kabi Container until it is totally empty. You may remove the air-vent needle for disposal.

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3 mm 3 mm

UK

Fresenius Kabi irrigation transfer sets are compatible with the Kabi Container: Bladder Syringes Filling Set (Ref 2841001), Transfer Set 1D (Ref 2883101), Transfer Set 2 D/X (2877201), Twin Transfer Set 7 D (2890501), AT1 Auto-transfusion Set (9005101), Cystoscopy Transfer Set (Ref 2883131).

Keep out of the sight and reach of children!
Do not use NaCl 0.9% Fresenius after the expiry date.

Date of revision
November 2021

Presentation
Polypropylene container (Kabi Container)
1 x 10 liters

Additional Information
Not for intravenous use or injection!
The sterility of NaCl 0.9% Fresenius is guaranteed only until the container is opened. The solution must be used immediately after opening.
Use only if the solution is clear and free from particles and if the package is undamaged.
For single use only.
Any unused portion of the solution must be discarded.
Re-use of medical single use equipment involves a risk of contamination of the equipment. This may cause infection in the patient or the user, which may lead to illness or death.

Seite 3

3 mm 3 mm

UK **INSTRUCTION FOR USE - READ CAREFULLY!** **Fresenius Kabi AG**
61346 Bad Homburg
Germany

NaCl 0.9% Fresenius - Solution for Irrigation

Composition
1 liter contains:
Sodium chloride 9 g
Hydrochloric acid (pH adjustment)
Sodium hydroxide (pH adjustment)
Water for injections in bulk, Ph. Eur.
Equivalent to: Na⁺ 154 mmol/l
Cl⁻ 154 mmol/l
pH value: 5.0 - 7.0
Titrations acidity < 1 mmol NaOH/l
Theoretical osmolality 308 mosm/l

Contraindications
None known.

Precautions
Use with caution in patients with hyperhydration, hyponatremia, and hyperchloremia.

Side-effects
Hyperhydration, hyponatremia, hyperchloremia and metabolic acidosis may occur in case of absorption of large amounts of NaCl 0.9% Fresenius.

Interaction with other drugs
None known.

Dosage
Dosage and quantity to be applied will result from the type of intervention.

Descriptions and connectivity
The Kabi Container is a ten-liter canister. Cap and stopper can either be completely removed for pouring the solution or an irrigation transfer set can be spiked into the stopper of the Kabi Container. Based on tests conducted the following

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Seite 6

3 mm 3 mm

Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

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Indicates the need for the user to consult the instructions for use.

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Indicates the manufacturer's batch code so that the batch or lot can be identified.

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Indicates the date after which the medical device is not to be used.

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Indicates the manufacturer's catalogue number so that the medical device can be identified.

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Seite 9

3 mm 3 mm

CE marking of conformity with Medical Devices Directive 93/42/EEC (as amended) with the intervention of Notified Body 0123.

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Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile.

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Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

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Indicates the medical device manufacturer.

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Check product name, container integrity, solution composition, lot number and product expiry date. Inspect the container for damage or leakage and for particles. Do not use if the container is damaged or if the solution contains particles. Do not use if tamper-evident plug has been removed or screw-seal is broken.

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Identify the plastic screw cap and the red rubber stopper underneath.

Tear off the tamper-evident plug/cover from the cap. Do not touch the surface of the rubber stopper!

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Seite 11

3 mm 3 mm

Indicates a medical device that should not be used if the package has been damaged or opened.

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This product is not for intravenous use or injection.

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Indicates a medical device that needs protection from light sources.

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Indicates a medical device that needs to be protected from moisture.

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Indicates the temperature limits to which the medical device can be safely exposed.

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