MEDICAL DEVICES

BioR Blood Filters

LEUKOCYTE DEPLETION FILTERS FOR RED CELL CONCENTRATES





BioR Blood Filters

Leukocyte depletion filters for red cell concentrates

BioR product range guarantees high filtration efficiency and performance, excellent recovery, user friendliness and time saving.

They are available as blood bank or bedside use filters in a wide variety of different configurations.

BioR 01 plus / max


Filter made by non woven material. Rigid, transparent housing, with printed batch number. Gamma rays sterilization. 3 years shelf life from date of manufacture.

Specifications

- Filtration capacity: 1 unit of red cell concentrates
- Filtration efficiency:
 - consistently averaging less than 0,1 x 10⁶ residual leukocytes (01 plus)
 - consistently averaging less than
 0,5 x 10⁵ residual leukocytes (01 max)
- Red blood cell recovery: averaging 88%*
- Filter housing hold-up volume: 33 ml*
- Filtration time: consistently less than 10 min (blood bank use)
- No priming with saline required
- Draining valves with integrated air filters (auto-venting) allow additional emptying of the filter housing and tubing to maximise the recovery of hemoglobin
- 200 µm prefilter for microaggregate removal
- Latex free

Leukocyte filters are characterized by:

- High filtration efficiency and performance
- Easy handling
- Minimal volume loss
- Short filtration time





* after use of auto

BioR 01 plus (n = 215) ²	PRE Volume (ml)	POST Volume (ml)	Htc (%)	Hb Recovery (%)	PRE WBC (x 10 ⁸ /unit)	POST WBC (x 10 ⁴ /unit)	Filtration Time (mm:ss)
Mean	282	250	61,1	88	4,0	3,5	08:12
Stand. Dev.	20,2	20,1	3,6	2,1	4,19	2,20	01:58
Min	203	175	50,0	83	0,3	1,8	04:26
Max	344	314	69,0	96	26,3	15,3	14:17

² Data from Center 1 (Italy, 2010): BC-depleted RCC units, stored in SAG-M at +4±2°C for 2-40 days; residual WBCs counting by Nageotte chamber.



The computer

provides optimal flow characteristics.

enables backtracking and follow-up

quality control. For optimal priming, wetting and filling of the filter the

The batch number printed on each filter

configuration is assembled in an upside down position. In case of emergencies blood can be filtered with a maximum pressure of 300 mm Hg or 0.4 bar.

Unique features of BioR

Filtering material¹

- Melt-blown non-woven polyester fibre
- Non ionic coated fibre surface
- Neither positively nor negatively charged
- Excellent wetting characteristic
- High biocompatibility

¹Trans. Med. Rev. 1995, 9 (2), 145-166

BioR 02 plus

Filter made by non woven material. Rigid, transparent housing, with printed batch number. Gamma rays sterilization. 3 years shelf life from date of manufacture.

Specifications

- Filtration capacity: 2 units of red cell concentrates
- Filtration efficiency: consistently averaging less than 0,2 x 10⁶ residual leukocytes
- Red blood cell recovery: averaging > 90% *
- Filter housing hold-up volume: 50 ml*
- Filtration time: consistently less than 20 min (for both units, blood bank use)
- No priming with saline required
- Draining valves with integrated air filters (auto-venting) allow additional emptying of the filter housing and tubing to maximise the recovery of hemoglobin
- 200 µm prefilter for microaggregate removal
- Latex free




-venting valves

BioR 01 max (n = 130) ²	PRE Volume (ml)	POST Volume (ml)	Htc (%)	Hb Recovery (%)	PRE WBC (x 10 ⁸ /unit)	POST WBC (x 10 ⁴ /unit)	Filtration Time (mm:ss)
Mean	281	249	60,7	88	3,5	3,4	08:35
Stand. Dev.	18,9	19,3	3,4	2,2	3,67	2,55	01:59
Min	233	201	50,1	83	0,2	2,0	05:10
Max	356	317	69,2	92	22,8	18,5	14:50

² Data from Center 1 (Italy, 2010): BC-depleted RCC units, stored in SAG-M at +4±2°C for 2-40 days; residual WBCs counting by Nageotte chamber.

Ordering information

For more information such as technical details and manuals, please contact your local sales representative or us.



A variety of other configurations is available.



Fresenius Kabi Deutschland GmbH Else-Kröner-Straße 1 · 61352 Bad Homburg · Germany Phone: +49 (0) 6172 686-4912 · Fax: +49 (0) 6172 686-4728 biorbiop@fresenius-kabi.com · www.fresenius-kabi.com Your local sales organization: