Document No.: 472312969
Revision: A
Change Order: C0-034051

Title: PRODUCT INSERT AMICUS APHERESIS KIT -

SINGLE NEEDLE K6R2312

Page: 1 of 4

Last Modified: 21 DEC 2018 MKF

Color: Black

PROOFREADING INSPECTION / RELEASED ARTWORK

Proofreading Approval:

Print Name Signature Date

Code REF K6R2312

3x

FRESENIUS

(en) AMICUS Apheresis Kit - Single Needle

For Extended Platelet Storage with Plasma and Red Blood Cell Collection. For Use with the AMICUS Separator.

Product should be visually inspected prior to use

Store at Controlled Room Temperature. Protect from freezing. Avoid excessive heat.

Definition of "Controlled Room Temperature":

"A temperature maintained thermostatically that encompasses the usual and customary working environment of 20° to 25°C (68° to 77°F); that results in a mean kinetic temperature calculated to be not more than 25°C; and that allows for excursions between 15°C and 30°C (59° and 86°F) that are experienced in pharmacies, hospitals, and warehouses. Provided the mean kinetic temperature remains in the allowed range, transient spikes up to 40°C are permitted as long as they do not exceed 24 hours ... The mean kinetic temperature is a calculated value that may be used as an isothermal storage temperature that simulates the non isothermal effects of storage temperature variations."

Reference: United States Pharmacopeia, General Notices. United States Pharmacopeial Convention, Inc. 12601 Twinbrook Parkway, Rockville, MD.

Description:

This apheresis kit consists of:

- (1) 1000 mL 0.9% Sodium Chloride Injection, USP (saline)
- (1) 1000 mL Anticoagulant Citrate Dextrose Solution, USP, Formula A (ACD-A)
- (2) 1000 mL plastic platelet storage containers
- (1) 800 mL plastic plasma storage container
- (1) 600 mL plastic red blood cell container
- (1) 17 gauge needle

Indications and Usage:

This kit is designed for use with the AMICUS separator for the collection of plasma, red blood cells, and leukoreduced platelets.

Platelet concentrate may be stored in the platelet storage container(s) as follows:

• For up to 5 days at 20°C to 24°C with continuous gentle agitation.

Each container can store up to a maximum of 4.7 x 10¹¹ platelets in the appropriate plasma volume. Refer to the Operator's Manual for minimum plasma storage volumes.

Best results will be achieved by measuring the donor's platelet count immediately prior to collection and entering this value into the Estimator.

The minimum post-procedure platelet count should be greater than or equal to $100,\!000/\mu L.$

We recommend not processing more than 5500 mL of whole blood for a single platelet product, 7000 mL of whole blood for a double platelet product, or 8000 mL of whole blood for a triple platelet product.

Plasma and red blood cells must be handled (freezing, storage conditions, labeling, etc.) according to the applicable procedures and regulations. Red blood cells must further be processed using a cleared AMICUS red blood cell set. Red blood cells should be stored in an appropriate storage solution and container.

Total plasma collected should not exceed 600 mL.

Plasma collected for use as fresh frozen plasma must be processed and frozen within 8 hours of successful venipuncture.

Contraindication:

Use of the AMICUS separator is contraindicated in those cases where adequate anticoagulation cannot be achieved.

Refer to the AMICUS Operator's Manual for a complete listing of contraindications, adverse effects, warnings and cautions.

Precautions:

- Procedures involving extracorporeal circulation are associated with the
 risks of blood loss, hemolysis, air embolism, and/or blood clotting. To
 minimize these risks, do not deviate from the procedures described in the
 Operator's Manual.
- To protect system integrity, do not prime the inlet/return line until immediately before venipuncture.
- 3. The separator should be operated only by a qualified healthcare professional who has been trained in the operation of this instrument.
- Failure to follow instructions prior to venipuncture may result in dilution of whole blood sample.
- The donor should meet appropriate regulatory requirements for the intended procedure and be in good health on the day of the donation.
- 6. Donors with low platelet counts may provide less than expected product yield and may have less than the usual post-donation platelet count. The platelet count of the donor should be taken into consideration when performing any platelet apheresis procedure. Because of increased platelet collection efficiency of the AMICUS separator, it may be necessary to reduce the volume of whole blood processed, compared to other blood cell separators, to avoid unacceptable low post-donation platelet counts.
- 7. Modification of the closed system or replacement of the needle under open system conditions negates any claims to a closed system, unless approved sterile connection methods are employed. If closed system conditions have been maintained, the platelet product may be stored for up to 5 days.

If open system collection conditions have occurred, platelets must be stored no longer than 24 hours. Any loss of sterile closure following collection results in a platelet product which must be stored no longer than 4 hours.

- 8. If the INTERLINK injection site on the inlet/return line is used, the closed system has been opened and platelets must be stored no longer than
- The use of a sterile connecting device to modify the configuration of the apheresis kit should be done using an operating procedure that is consistent with applicable regulatory guidelines and the manufacturer's directions for use.
- 10. Refer to the Operator's Manual for conditions when platelet products should be counted for residual white blood cells.

-1- 47-23-12-969 REV: A

Document No.: 472312969 Revision:

Change Order: CO-034051

PRODUCT INSERT AMICUS APHERESIS KIT -Title:

SINGLE NEEDLE K6R2312

Page:

Last Modified: 21 DEC 2018 MKF

Color: Black

PROOFREADING INSPECTION / RELEASED ARTWORK

Proofreading Approval:

Print Name Signature Date

Adverse Effects:

Adverse donor reactions, identical to those seen with routine blood collection procedures, may occur. These include "lightheadedness," fainting, vomiting, hyperventilation, and hematoma formation at the site of phlebotomy. Syncopal reactions due to hypovolemia may occur.

Reactions unique to apheresis collection procedures may also occur. Allergic symptoms including skin redness, itching, hives, etc., may be seen. Cooling of saline or donor blood may be sufficient to induce chills. Infusion of anticoagulants containing citrate may lead to donor symptoms of moderate hypocalcemia due to chelation of calcium by unmetabolized citrate. Such donor reactions are usually manifested by a "tingling" feeling, often around the mouth or fingers. Other manifestations may include muscle discomfort, muscle twitching or spasms, and/or the presence of an unusual taste in the mouth. Return of the donor's blood should be temporarily slowed by reducing the citrate infusion rate (CIR) or the procedure may be halted if symptoms occur.

Symptoms of severe hypocalcemia rarely seen in donors undergoing leukapheresis include tetany, seizure, cardiac arrhythmia, and death.

Complications such as blood loss, hemolysis, air embolism, and/or blood clotting, may be associated with improper operating conditions.

This container is of PVC material that used di-(2-ethylhexyl) phthalate (DEHP) as plasticizer, where DEHP is reported to affect the development of male sexual organs and sperm formation in young animal tests. Although infinitesimal amount of DEHP may be eluted in such PVC containers, there is very little or no risk from exposure to DEHP. Hence, as the risk that may appear from not using this medical device is much larger than the risk concerned due to DEHP, there is no need to avoid using it.

Directions for Use:

Refer to Operator's Manual for the AMICUS separator for the installation procedure and complete instructions for use and plasma volumes for platelet

Dispose of waste in appropriate biohazard container or according to local regulatory requirements.

Symbols with Definitions:

Caution, consult instructions for use

Sterilized by a combination of steam and radiation. Sterile fluid path.

Non-pyrogenic fluid path Do not reuse

Do not use if the product sterile barrier system is compromised

Do not vent

Contains or presence of phthalate: Diethylhexylphthalate (DEHP)


LOT

Code Lot

Expiration date

Tare Weights:

Below are the approximate tare weights of the platelet container and associated kit components. Refer to the appropriate illustration for the tare weights for platelet product calculations.

	Approximate Weights (grams)	
Two (2) platelet storage containers with: Two (2) sample pouch assemblies, Two (2) roller clamps and tubings, One (1) Y-Junction, and One (1) 41.9 cm tubing	77	
Two (2) platelet storage containers with: Two (2) sample pouch assemblies, Two (2) roller clamps and tubings, and One (1) Y-Junction	73	
One (1) platelet storage container with: One (1) sample pouch assembly, One (1) roller clamp and tubing, and One (1) Y-Junction	37	
One (1) platelet storage container with: One (1) sample pouch assembly	31	
One (1) platelet storage container	27	
Plastic tubing	0.12 g/cm	
One (1) plasma storage container (800 mL)	39	
One (1) red cell container (600 mL)	32	

Fresenius Kabi AG

61346 Bad Homburg / Germany Tel.: +49 (0) 61 72 / 686-0 www.fresenius-kabi.com

Fenwal International. Inc.

Carretera Sanchez Km 18.5 Parque Industrial Itabo Zona Franca Ind. De S.C. Haina, Dominican Republic

INTERLINK is a trademark of Baxter International, Inc. © 2018 Fresenius Kabi AG. All rights reserved.

-2-

Document No.: 472312969
Revision: A
Change Order: C0-034051

Change Order: CO-034051
Title: PRODUCT INSERT AMICU

PRODUCT INSERT AMICUS APHERESIS KIT -

SINGLE NEEDLE K6R2312

Page: 3 of 4

Last Modified: 21 DEC 2018 MKF Color: Black **Proofreading Approval:**

Print Name Signature Date

PROOFREADING INSPECTION / RELEASED ARTWORK

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