# Ambix nova

## Inspirert av pasienters liv

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### Infusion

#### Infusion modes

- Flow rate only mode: continuous infusion at a flow rate
- Continuous mode: continuous infusion at a flow rate, volume and time
- Ramp mode: 3 phases infusion at a volume and time
- Intermittent mode: cyclic infusion at a flow rate, volume and time

#### Flow rate

From 1 mL/h to 600 mL/h, 1 mL/h increments

#### Accuracy

+/- 5% for flow rate +/- 10 mL/h +/- 0.5 mL for flow rate < 10 mL/h

#### Volume range

From 1 mL to 9999 mL, 1 mL increments

#### Priming

Automatic, semi-automatic and manual priming are available

#### Delayed start

Delayed start enables to postpone start of a new infusion from 1 min. to 24 h

#### KVO (keep vein open)

KVO function is a low flow rate from 1 mL/h to 10 mL/h to keep vein open during all infusion interruptions

#### Night mode

Night mode enables to manage pump luminosity at night by decreasing screen brightness

#### Max flow rate

Customization of pump flow rate to not exceed (from 1 mL/h to 600 mL/h)

#### Keypad lock

Disabling keypad to prevent unintentional key press

#### Nurse setting lock

Possibility to prevent patient from changing settings, access secured by code

### Alarms / Informations

#### High priority alarms (Pumping stops and red LED blinking)

Missing set, Door open, Upstream occlusion, Downstream occlusion, Empty bag / Air inline, End of infusion, Empty battery, Technical alarms

#### Low priority alarms (Pumping continues and orange LED is fixed)

End of infusion pre-alarm, Empty battery pre-alarm, Start reminder

#### Occlusion alarm

2 levels:

Low: 487.5 mmHg +/- 187.5 mmHg

High: 787.5mmHg +/- 262 mmHg

#### Air alarms

Single bubble volume setting: from 0.125 mL to 0.75 mL

Cumulated bubbles volume setting: from 0.25 mL to 0.75 mL

Duration of cumulated bubbles detection: from 1 min to 30 min.

#### End of infusion pre-alarm

Customizable alarm to notify approaching end of infusion from 1 min. to 59 min. or from 1 mL to 999 mL

#### History menu

Infusion history (250 infusion events), Cumulative counter of volume infused during a period of time (from 1 to 99 days), Alarm history (250 alarm events)

#### Reset manufacturing settings

Automatic reset of all settings and history information

### Technical specifications

#### Pumping mechanism

Linear peristaltic pumping system

#### Display

Color LCD 64.8 mm x 48.7 mm (320 pixels x 240 pixels)

#### Dimensions (H/W/D) / Weight

Pump: Approx. 156.4 mm x 129.7 mm x 49 mm / ~ 0.66 kg

Smart Holder Power: Approx. 132 mm x 118 mm x 46 mm (without pole clamp) / ~ 0.45 kg

Ambix is a registered trademark by Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most updated information. These medical devices are regulated by health authorities. Please refer to the respective CE Declaration of Conformity for the applicable Medical Device Directive or Regulation. Pictures are not contractual.

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