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Stand der Information

August 2012

PACKAGE INSERT – READ CAREFULLY!**Fresenius Kabi AG**

Medical Devices
D-61346 Bad Homburg v. d. H.

GB**Isotonic Saline**

Irrigation Solution

Composition	1 litre contains:
	Sodium chloride 9,0 g
	Na ⁺ 154.0 mmol/l
	Cl ⁻ 154.0 mmol/l
	Aqua inieciabilia
	Theoretical osmolarity 308 mosm/l

Indications Plasmaitonic (physiological) saline for internal and external application such as perfusion of the extracorporeal system in haemodialysis, postoperative bladder irrigation in all types of urological surgery, irrigation of the gastro-intestinal tract, fistulas and drainages. Also in wound treatment and to moisten cloths and dressings.

Contraindications Hypernatraemia, hyperhydration.
Note: Caution in case of hyperchloraemia.

Side-effects Not known when administered as directed.

Interaction with other drug Not known.

Application and Dosage

If not otherwise prescribed the dosage depends on its internal or external use according to the surgical intervention.

Presentation

Plastic-bags	4 x 3000 ml
Plastic-bags	2 x 5000 ml

Additional Information

Plastic containers are sometimes damaged on their transport from the manufacturer to the hospital or within the hospital. This may cause a contamination of the solution with microbial or fungal growth.

It is therefore recommended that a thorough visual inspection is carried out prior to connecting the container and using the solution. Attention must be paid to the seal, seam and corners of the container.

Isotonic Saline must not be used after expiry date!

Use only clear solutions and undamaged containers.

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