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**Stand der
Information**

August 2012

PACKAGE INSERT – READ CAREFULLY!**Fresenius Kabi AG**

Medical Devices
D-61346 Bad Homburg v. d. H.

**Purisole® SM diluted**

Irrigation Solution

Composition

1 litre contains:

Sorbitol	27.0 g
Mannitol	5.4 g
Water	
pH-value	4.5 – 7.0
titration acidity	< 1 mmol NaOH/l
theor. osmolarity	178 mosm/l

Indications

Intra- and postoperative bladder drainage in urological surgery such as transurethral resection of the prostate or bladder tumors, suprapubic examinations, cystoscopic examinations, lithotripsy, among others.

Distension and irrigation of the uterine cavity in operative hysteroscopy, e.g. endometrial ablation, septum dissection, resection of submucous myoma and synechiolysis.

In diagnostic procedures, e.g. diagnostic hysteroscopy and cystoscopy.

Contraindications

Methanol intoxication. Infants and small children. Known hypersensitivity to Sorbitol and / or Mannitol.

Purisole® SM diluted contains sorbitol. For this reason, the presence of a hereditary fructose intolerance or a lack of fructose 1–6 diphosphatase must be ruled out prior to irrigation with Purisole® SM diluted. Should anamnestic data be unavailable or unclear, a fructose intolerance test must be performed before irrigation. Non-observance of these contraindications may lead to potentially fatal complications.

Side effects A possible intravascular absorption of large amounts of sodium-free irrigation solution may result in changes in electrolyte balance (hyponatremia) and, as a consequence, to a change in acid base balance (acidosis).

Agitation, nausea, dyspnea, as a sign of hyperhydration of the patient, are early symptoms of a imbalance (TUR syndrome).

Interactions with other medicines Not observed.

Dosage instructions Dosage and quantity to be applied result from the type of intervention.

Method of application Do not use for infusion or injection!
For irrigation.

Presentation	Plastic-bags	6 x 1500 ml
	Plastic-bags	4 x 3000 ml
	Plastic-bags	2 x 5000 ml

Additional Information Plastic containers are sometimes damaged on their transport from the manufacturer to the hospital or within the hospital. This may cause a contamination of the solution with microbial or fungal growth.

It is therefore recommended that a thorough visual inspection is carried out prior to connecting the container and using the solution. Attention must be paid to the seal, seam and corners of the container.

In case of doubt, the attending physician should decide on the use of the solution.

Use the solution only if container and seal are undamaged and the solution has been released by the attending physician.

Use only clear solutions and undamaged containers.

Purisolet® SM diluted must not be used after expiry date!

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