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Stand der Information

September 2012

PACKAGE INSERT – READ CAREFULLY!**Fresenius Kabi AG**

Medical Devices
D-61346 Bad Homburg v. d. H.

Ringer's Solution for Irrigation

Solution for irrigation and wound lubrication

Composition

1 litre contains:

sodium chloride	8,6 g
Potassium chloride	0,3 g
Calcium chloride 2 H ₂ O	0,33 g
Na ⁺	147,2 mmol/l
K ⁺	4 mmol/l
Ca ⁺⁺	2,25 mmol/l
Cl ⁻	155,7 mmol/l
Water	
pH-value	5,0–7,0
titrations acidity	< 1 mmol NaOH/l
Theoretical osmolarity =	309 mosm/l

Indications

Irrigation solution for internal and external use:
During operative procedures to rinse the area being operated on and to keep body tissues moist.

For the lubrication of surgical tampons and dressings.

For the irrigation of wounds in cases of external trauma and burns.

For irrigation during diagnostic examinations.

Contraindications

Hypernatraemia, hyperhydration.

Note: Caution in case of hyperchloraemia.

Side-effects

None have been reported.

Interaction

None have been reported.

Dosage and Application

For internal and external irrigation.
For the lubrication of wounds and surgical dressings.

Unless otherwise prescribed, the amount of solution required corresponds to the area of application.

Note:

Do not use for intravenous injection or infusion. The sterility and apyrogenicity of Ringer's Solution are guaranteed only until the opening of the container. The contents must be used at once. Do not keep any remaining solution for later use. Discard all left-over solution.

Presentation

Plastic-bags	6 x 1500 ml
Plastic-bags	4 x 3000 ml
Plastic-bags	2 x 5000 ml

Additional Information

Plastic containers are sometimes damaged on their transport from the manufacturer to the hospital or within the hospital. This may cause a contamination of the solution with microbial or fungal growth.

It is therefore recommended that a thorough visual inspection is carried out prior to connecting the container and using the solution. Attention must be paid to the seal, seam and corners of the container.

Keep out of the reach of children.

Do not use Ringer's Solution for Irrigation after the expiry date.

Use only if the solution is clear and the container undamaged.

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