



MEDICAL DEVICES
Innovating to ease the flow of life

Disposables Information



Syringe Pump
Sets


Name: Injectomat Line ASV Product code: 9004149



- Syringe pump line
- ⇒ Line with integrated free-flow protection for syringe pumps
 - ⇒ With Anti-Siphon valve on female luer port
 - ⇒ Additional Robson clamp for additional security during longer infusion stops or when the syringe itself is manipulated.
 - ⇒ Contains DEHP



Latex free



Specifications:

- Tube Material:
PVC ID 0.95mm
- Length:
200 cm
- Priming volume:
approx 1.4 ml
- Anti-Siphon Valve:
opening pressure
3psi ±1.5psi
0.2bar ±0.1bar
- Units:
200 per box
- Pressure resistance:
2 bar




MEDICAL DEVICES
Innovating to ease the flow of life

Disposables Information



Syringe Pump
Sets

Name: Injectomat Line TIVA Twin Product code: 9004155



Dual syringe pump line for the use in intravenous anesthesia

- ⇒ Dual line for the easy and safe application of both drugs used in TIVA anesthesia
- ⇒ With Back-Check valve on all female luer ports to direct the drug flow
- ⇒ Additional Robson clamp for additional security during longer infusion stops or when the syringe itself is manipulated.
- ⇒ Contains DEHP



Latex free



Specifications:

- Tube Material:
PVC ID 0.95mm
- Length:
210 cm
- Priming volume:
approx 2 x 1.4 ml
+ 0.04ml
- Back-Check Valve:
opening pressure
0,27 psi / 0.018 bar
- Units:
200 per box
- Pressure resistance:
2 bar




MEDICAL DEVICES
Innovating to ease the flow of life

Disposables Information



Syringe Pump
Sets

Name: Injectomat Line TIVA Twin Y Product code: 9004163



Dual syringe pump line with additional infusion port for the use in intravenous anesthesia

- ⇒ Dual line for the easy and safe application of both drugs used in TIVA anesthesia
- ⇒ With Back-Check valve on all female luer ports to direct the drug flow
- ⇒ Additional Robson clamp for additional security during longer infusion stops or when the syringe itself is manipulated.
- ⇒ Contains DEHP



Latex free



Specifications:

- Tube Material:
PVC ID 0.95mm
- Length:
210 cm
- Priming volume:
approx 2 x 1.4 ml
+ 0.08ml
- Back-Check Valve:
opening pressure
0,27 psi / 0.018 bar
- Units:
200 per box
- Pressure resistance:
2 bar