



**FRESENIUS
KABI**

caring for life

Freka[®] Belly Button

PERCUTANEOUS LOW PROFILE FEEDING TUBE
For Long-Term Enteral Feeding



Freka Belly Button

Percutaneous low profile feeding tube for long-term enteral feeding

Small optimised design




> For more discretion and patient comfort



- › FR/CH size and length are labelled for easy identification
- › Balloon fill volume printed on button head for user reassurance
- › Range of accessories available

Freka Belly Button pack contains:

- › Freka Belly Button
- › 1 Introducer
- › 1 Syringe for balloon filling
- › 10ml Ampoule sterile water
- › 1 Lubricant gel and 2 gauze pads

Freka Extension Set 30cm	Introducer	Freka Stoma Measuring Device
Product code: 7981397	To support the shaft for ease of placement*	Product code: 7750921
 A clear plastic extension set with a purple connector and a white syringe-like bulb.	 A long, thin, green plastic introducer with a bulbous head.	 A white plastic measuring device with a long shaft and a circular head.

* For detailed placement refer to instruction manual

Ordering Information



Freka Belly Button	FR size/ shaft length	FR 12	FR 14	FR 16	FR 18	FR 20	FR 24
Containing: 1 Freka Belly Button 1 Introducer (simple) 1 Syringe for balloon filling 10ml Ampoule sterile water 1 Lubricant gel 2 Gauze pads	1.0 cm	7751612	7751622	7751632	7751642	7751652	7751662
	1.5 cm	7751613	7751623	7751633	7751643	7751653	7751663
	2.0 cm	7751614	7751624	7751634	7751644	7751654	7751664
	2.5 cm	7751615	7751625	7751635	7751645	7751655	7751665
	3.0 cm	7751616	7751626	7751636	7751646	7751656	7751666
	3.5 cm	7751617	7751627	7751637	7751647	7751657	7751667
	4.0 cm	7751618	7751628	7751638	7751648	7751658	7751668
	4.5 cm	7751619	7751629	7751639	7751649	7751659	7751669
	5.0 cm	7751620	7751630	7751640	7751650	7751660	7751670

Accessories

	Art.No.
Freka Button Extension Set 1, 30 cm	7981397
Freka Stoma Measuring Device	7750921