



AmiCORE

Apheresis System

A new level of
platelet collection

Easy to operate, ergonomic design

Smart monitoring technology

Enhanced donor management tools

Thoughtful design and exceptional
for a high level of efficiency

SIMPLE

AmiCORE's innovative design
streamlines routine operations

- **SmartRoute** simplifies kit installation with a tray holder, a single cassette, tubing guides and color coded tubing
- **Centrifuge** positioned at optimal height to streamline kit installation
- **Intuitive User Interface** provides easy navigation with large, user-friendly touch screens
- **Scales** provide accurate fluid volume measurements
- **Single Needle** procedure for donor comfort

AmiCORE helps to simplify each step of the platelet collection process, allowing operators more time to focus on donors.



performance come together
and donor management.

SMART

AmiCORE's smart monitoring brings your technology up to date

- **TriVision** provides precise optical monitoring of the cellular interface position to enhance separation
- **Platelet Post-Count Accuracy** allows selection of target yields consistent with donor safety
- **Process Leukoreduction** automatically provides leukoreduced platelets without the additional step of filtration
- **Electronic Procedure Records** simplify the implementation of GMP documentation

AmiCORE provides an extensive set of donor management tools that operators can rely on

- **Saline Infusion** replaces fluid during donation
- **Intelligent Flow Control (IFC)** automatically establishes optimal donor flow rates
- **Real-time Pressure Monitoring and Numerical Flow Rate Control** help protect the donor's vein and maintain comfort
- **Built-in Pressure Cuff** prompts donors to squeeze during the draw cycle and relax during the return cycle



AmiCORE Apheresis System

Simple and Smart

AmiCORE's enhanced donor management tools provide blood centers with the flexibility to customize the platelet collection procedure to target a wide range of donors.



Ordering Information

6R8800	AmiCORE Apheresis System
X6R8880 P6R8880	AmiCORE Apheresis Kit, single needle
X6R8882 P6R8882	AmiCORE Apheresis Kit, single needle with two platelet containers
6212836539	AmiCORE USB Flash Drive

Refer to AmiCORE Operator's Manual for a full list of warnings and cautions associated with the use of the AmiCORE device.

SmartRoute and TriVision are trademarks of Fresenius Kabi AG.