




# Overview of ENFit Enteral Pump Sets

**Amika<sup>+</sup>**



	EasyBag		VarioLine		Bag		Bag Mobile	
--	---------	--	-----------	--	-----	--	------------	--

ART. NR.	7751907	7751900	7751909	7751919	7751715T	7751914	7751906	7751915
<b>UPPER END</b>								
ENPlus spike	●	●	●	●				
VarioLine bottle adapter			●	●				
Attached Bag					●	●	●	●
<b>LINE COMPONENTS</b>								
Drip chamber	●	●	●	●	●	●		
Roller clamp							●	●
ENFit medication port	●	●	●	●	●	●	●	●
Cover on medication port				●		●		●
<b>LOWER END</b>								
Adaptors for ENLock and LuerLock	●		●					
ENFit connector	●	●	●		●		●	
ENFit connector with cover				●		●		●
Tubing length (mm)	250	250	250	250	250	250	242	242
Length med. port to lower end (cm)	40	40	40	40	40	40	40	40

Not made with DEHP

\* Pump set with transparent cover on the lower-end ENFit connector and tethered cover on the ENFit medication port