




**FRESENIUS
KABI**

caring for life

Freka[®] Jejunostomisæt med ENFit



Freka Jejunostomisæt:

Polyuretansonde

Røntgenstribe og interne
sondemarkeringer hver 10. cm

Foldbar silikone fikseringsplade

CH 9 / 75 cm

Varenr. 7755645

Freka® Jejunostomisæt med ENFit

Article Freka® FCJ FR 9, ENFit
Article No. FK 7755645 (M90800308)

Parts List

Pos.	Description	Material
1	Tube 1.85 x 0.55 x 750 ± 5 mm with one x-ray contrast strip	PUR (Polyurethane) and PUR + 60% BaSO ₄
2	Robert's clamp, white	POM (Polyoxymethylene)
3	External fixation plate	Polydimethylsiloxane (Silicone elastomere), Bariumsulphate
4	Fixation clam, white	PBT (Polybutylenterephthalate), Colourant
5	ENFit male - PEG Connector, white	PBT (Polybutylenterephthalate), Colourant
6	Strap, transparent	LDPE (low density polyethylene)
7	ENFit female cap, lilac	ABS (Acrylonitrine-Butadiene-Styrene), Colourant
8	Mini tube 2 x 0.18 x 11.5 mm	Stainless steel
9	Fixation ring	SEBS (Styrene-Ethylene/Butylene-Styrene)
10	Fixation screw	PBT (Polybutylenterephthalate), Colourant
11	Fixation aid	PBT (Polybutylenterephthalate), Colourant
12	Partable cannula 3.55 x 80 mm, blue	Protective cap: PP (Polypropylene) Hub for cannula, moulded: ABS (Acrylonitrine-Butadiene-Styrene), Colourant Cannula: Stainless steel Silicone oil (Polydimethylsiloxane)
13	Plug for cannula mandrin, white	MABS (Methylmethacrylate-Acrylonitrile-Butadiene-Styrene), Colourant
14	Partable cannula 3.55 x 80 mm, white	Protective cap: PP (Polypropylene) Hub for cannula, moulded: ABS (Acrylonitrine-Butadiene-Styrene), Colourant Cannula: Stainless steel Silicone oil (Polydimethylsiloxane)
15	Cannula mandrin 3 x 154 mm	Stainless steel

EN_0073/Maj 2017