



**FRESENIUS  
KABI**

caring for life

# Freka<sup>®</sup> Tragtadapter med ENFit



**Freka Tragtadapter:**

---

**med ENFit**

---

**til dekompression med sug**

---

**for at skabe kompatibilitet  
mellem sonder med ENFit og  
enteralsæt og spørøjter med  
ENLock**

---

**varenr. 7755695**

---

# Freka® Tragtadapter med ENFit

**Article** Freka® ENFit Universal Funnel Adaptor, ENFit  
**Article No. FK** 7755695 (M90800302)

## Parts List

Pos.	Description	Material
1	Universal funnel adaptor with Luer neg.	PUR (Polyurethane), Colourant
	ENFit Female - LL Male Connector	ABS (Acrolonitrile-Butadiene-Styrene)
<b>Packaging</b>		
	Bag for 15 pcs	LD-PE (Low Density Polyethylene)
	Label sticker 148 x 99 mm	Paper
	Carton 290 x 187 x 205 mm for 600 pcs.	Corrugated cardboard
	Label sticker 148 x 99 mm	Paper

EN\_0071/Marts 2017