

Freka® EasyBag Bolusadapter med ENFit

Article EasyBag Bolus Adapter, ENFit
Article No. FK 7755694 (M90800301)

Parts List

| Pos. | Description | Material |
|----------------------|---|---|
| 1 | Protective cap for ENPlus | PP (Polypropylene) |
| 2 | ENPlus spike | ABS (Acrylonitrile-Butadiene-Styrene), Colourant |
| 3 | Retainer nut for ENPlus | PP (Polypropylene), Colourant |
| 4 | Sealing for cross spike | PE (Polyethylene) |
| 5 | Adapter EasyBag spike | PUR (Polyurethane) |
| 6 | NoDop-tube 3 x 5 x 50 ± 5 mm | PVC (Polyvinylchloride) with plasticizer TEHTM |
| 7 | ENFit male | ABS (Acrylonitrile-Butadiene-Styrene), Colourant |
| 8 | Strap for ENFit cap | PP (Polypropylene), Colourant |
| 9 | ENFit female cap | ABS (Acrylonitrile-Butadiene-Styrene), Colourant |
| | Solvent | THF (Tetrahydrofurane) |
| | Solvent | CHN (Cyclohexanone) |
| Sterilization | | |
| | ETO / CO ₂ mixture | Ethylene oxide / Carbon dioxide |
| Packaging | | |
| | Compound foil 70 µm | PA/PE (Polyamide / Polyethylene) |
| | Paper 148 x 240 mm | Paper |
| | Carton 249 x 148 x 150 mm, for 15 pcs. | Corrugated cardboard |
| | Label sticker 145 x 94 mm | Paper |
| | Indicator adhesive tape | Tape |

EN_0072/April 2017