

AmiCORE

A New Level of Platelet Collection

SmartRoute

SmartRoute simplifies kit installation with a tray holder, a single cassette, tubing guides and color coded tubing.




Tray holder places the kit in front of the operator for easy installation.



Intuitive user interface provides easy navigation with large, user-friendly touch screens.



Single cassette installation saves time and simplifies the process.



AmiCORE Quick Report



AmiCORE quick reports simplifies the implementation of GMP documentation

Procedure Summary Report

707295
07:41 22/12/2015

00:36 Proc. Time	06:50 - 07:27	2.5 x10 ¹¹ Platelet Yield	21 mL RBC Loss
1726 mL WB Processed	172 mL Storage Fluid	26 mL AC in Product	90 mL Plasma Loss
209 mL AC Used			
288 mL Saline Used			

AmiCORE Quick Report*

FRESENIUS KABI

AmiCORE Single Donor Platelet Collection Form

Hospital/Blood Center Name: Blood Center

Operator ID: 707295 Procedure Type: Single-Needle Platelets in Plasma

• Donor Information

Donor name :	Donation ID : 429851214	Date : 2015-12-22
Donor ID: 567479832	HCT or Hgb : 43%	Height : 173 cm
Blood Type:	Gender : M	Body weight : 78 kg

• Disposable Kit and solutions

Disposable Kit code	Disposable Kit Lot No.	Kit Expiration Date
6R8882	FA16823158	2018-02-27
Saline code	Saline Lot No.	Saline Expiration Date
2295121E	14IF7314	2018-04-29
ACD code	ACD Lot No	ACD Expiration Date
487898Q	C991109	2017-03-30

• Platelet collection

	Product Yield Xe11		Donor Count K/ μ L	
	Target Yield		Pre-Count	
Pre procedure		2.5		250
Post procedure	Estimated Yield	2.5	Est. Post-Count	209

Refer to AmiCORE Operator's Manual for a full list of warnings and cautions associated with the use of the AmiCORE device.

*This report depicts sample information only. The information shown is used to express information available in an AmiCORE quick report.

