




**FRESENIUS
KABI**

caring for life

agilia
Connect



Vigilant® Bridge

EMR Auto-documentation Software

An auto-documentation interface for
PDMS connectivity

Reduce nurse charting time

Provide interoperability with PDMS/EMR

Sustainable investment by using HL7
standard and IHE profiles

Vigilant® Bridge

Product characteristics

Features: allows auto-documentation with PDMS/EMR using:

- HL7 V2: 6 messages
- IHE Patient Care Device (PCD) Technical Framework
- With Device Enterprise Communication (DEC) profile
- With Infusion Pump Event Communication (IPEC) profile

Link+ Agilia rack communicates with Centerium over the hospital 802.3 based intranet.

Language: English.

Hardware requirements

Vigilant Bridge software must be installed on the same server as Centerium Data Management Solution.

Please refer to Centerium data sheet for Software requirements, adding for Vigilant Bridge:

Software requirements

Vigilant Bridge software must be installed on the same server as Centerium Data Management Solution.

- OpenVPN server version 2.3.6
- Visual C++ Redistributable package for Visual Studio 2012 Update 4, 64 bit
- Python programming language version 2.6.5

All these software: Redistributable Package are installed automatically during software installation.

Connected devices requirements

Vigilant Bridge software is compatible with the following devices:

- Link 4+ Agilia (4 pump slots) with firmware version 3.0
 - Link 6+ Agilia (6 pump slots) with firmware version 3.0
 - Link 8+ Agilia (8 pump slots) with firmware version 3.0
- And
- All Agilia Connect VP pumps with firmware version 1.4 and 2.2.c
 - All Agilia Connect SP pumps with firmware version 1.6, 1.6a and 2.2

Pumps must be connected through a Link+ Agilia. A maximum of 100 Link+ can be connected simultaneously.

- Centerium 1.0b and 1.0c.

Cable requirements

- Ethernet RJ45 cable. Category 5e minimum.

Security requirements

Vigilant Bridge software must be installed on the same server as Centerium Data Management Solution.

Please refer to Centerium Administrator's Guide for Security requirements.

Storage area requirements

Vigilant Bridge software must be installed on the same server as Centerium Data Management Solution.

Please refer to Centerium data sheet (CC11077 datasheet Centerium) for Hardware and Software requirements.

Agilia and Vigilant are registered trademarks by Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most updated information.

