



Vigilant[®] Insight

Infusion Data Reporting Software

Monitor drug library adoption
to support your safety goals

Improve Drug usage and
practice

Decrease human errors and
limit alarm fatigue

Latest cybersecurity standards

Vigilant® Insight

Vigilant Insight is an application hosted in Centerium with the following prerequisites:

Hardware Requirements

The hardware requirements for Centerium are dependent on the number of pumps used, as more powerful hardware is needed when Centerium is used with a large number of pumps. Centerium was validated against the following host computer requirements:

	1 - 1000 pumps	1000 - 3000 pumps
Operating System	Microsoft Windows Server 2008 R2 SP1 (Version 6.1) Microsoft Windows Server 2012 (Version 6.2)	Microsoft Windows Server 2008 R2 SP1 (Version 6.1) Microsoft Windows Server 2012 (Version 6.2)
Configuration	2-tier	2-tier

Application Server

	1 - 1000 pumps	1000 - 3000 pumps
Cores	16	16
Processor	18 GHz	18 GHz
Total Physical Memory	64 GB	64 GB
Disk Size	200 GB	200 GB
Network Interface Card	IEEE 802.3 10/100/1000Mbps card	IEEE 802.3 10/100/1000Mbps card

Database Server

	1 - 1000 pumps	1000 - 3000 pumps
Cores	16	16
Processor	18 GHz	18 GHz
Total Physical Memory	32 GB	32 GB
Disk Size	200 GB	1000 GB
Network Interface Card	IEEE 802.3 10/100/1000Mbps card	IEEE 802.3 10/100/1000Mbps card
	Windows Server 2012	Windows Server 2008 R2 SP1
SQL Server 2012 Standard	Validated	Validated
SQL Server 2008 R2 SP2 Standard	Not supported	Validated

Security requirements caring for life

Vigilant Insight software must be protected with an up-to-date virus scanner.

For any other OS version or Database version:

The installer will display a message stating that the version has not been validated. Use with other versions must be validated prior to use.

The Centerium Installer enforces usage for specified versions; use with other versions will be prevented by the installer.

Vigilant is registered trademark by Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most updated information.

