



# Vigilant<sup>®</sup> Sentinel

## Your smart infusion status

A clinical application for the remote and centralized viewing of infusion pump status

---

Optimize your workflow

---

Reduce alarm fatigue

---

Manage your time

---

# Vigilant<sup>®</sup> Sentinel

## Product characteristics

**Features:** aggregated view at the ward level, detailed views at the bedside level and empty container view.

**Scalability:** 24 beds, 2 links per bed, 384 pumps.

**Languages:** English, Spanish, German, French, Dutch, Italian, Czech, Slovak, Brazilian Portuguese, Bahasa, Polish, Hungarian, European Portuguese.

## Hardware requirements

### Computer requirements:

Desktop PC mandatory, no laptop allowed.

### Display requirements:

- Display with resolution 1920 x 1080.
- Minimum recommended size 24 inches.
- Possibility to use touchscreen or mouse.

### Cable requirements:

Ethernet RJ45 cable: category 5e minimum.

## Computer software requirements

- Microsoft Windows 8.1, 64 bits.
- Microsoft Windows 10, 64 bits.
- .Net framework version 4.6.2 and Microsoft® Visual C++ 2015 Redistributable Package are installed automatically during software installation.
- Adobe Acrobat Reader.

Please refer to software Installation Instructions for additional information.

## Connected devices requirements

Vigilant Sentinel software is compatible with the following devices:

- Volumat standard with firmware version equal to 1.2e
- Volumat MC with firmware version equal to 2.3f
- Injectomat standard with firmware version equal to 2.5g
- Injectomat MC and TIVA with firmware version equal to 4.3h
- All VP pumps with firmware version 1.4 and later versions up to and including 2.2.x
- All SP pumps with firmware version 1.6a and later versions up to and including 2.2.x
- Link+ 4 (4 pump slots) with firmware version 3.0
- Link+ 6 (6 pump slots) with firmware version 3.0
- Link+ 8 (8 pump slots) with firmware version 3.0

## Security requirements

- Vigilant Sentinel software must be protected with an up-to-date virus scanner.
- Vigilant Sentinel software must be configured to support user authentication as described in the software Installation Instructions.

## Storage area requirements

A storage area must be available and allocated with the following specifications: Data volume: 50 MB.

## Delivery Media

- CD-ROM including software and product documentation in electronic and printed form.
- Key 2 download site at <https://key2.fresenius-kabi.com>

Agilia and Vigilant are registered trademarks of Fresenius Kabi in selected countries. Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for the most up-to-date information.

