



CompoSeal Slim

Tube sealing system

CompoSeal Slim provides a flexible desktop tube sealing system for perfect and reliable seals.

High quality seals

Reliable operations

Quick and easy handling

Cascade possibility (segment sealing)

CompoSeal Slim

Flexible desktop tube sealing system for single and cascade sealing

The CompoSeal Slim has the ability of cascade sealing up to 8 devices

- Segment length adjustable between 60 - 120 mm
- Serial sealing from start of segment tube towards blood bag, indicated by a light indicator
- Per two machines, only one power cord needed





High quality
seals for reliable
operations

Efficient, quick
and easy handling

Ability to cascade
sealing of up
to 8 devices

CompoSeal Slim ensures seal quality

- Increased range of sealable tube diameter (up to 5.9 mm)*
- Integrated spill protector
- Detachable sealing head for easy cleaning and disinfection
- Sealing time 1 -2 seconds
- Indication by light of sealing process and sealing errors
- Additional option for hand sealer

*Not applicable for the optional hand sealer

Optional hand sealer function

With an optimal ergonomic design the hand sealer requires very low hand force; visual and acoustic alarm; detachable head is easy to clean.

*CompoSeal Slim hand sealer is
available on request*



CompoSeal Slim Tube sealing system



Process efficiency to help you achieve more

Ordering information

For more information such as literature, technical details, working instructions and equipment, please contact your local sales representative.

CE marked, Class I according Directive 93/42/EEC, as amended by 2007/47/EC.



**FRESENIUS
KABI**

caring for life

Fresenius Kabi AG
61346 Bad Homburg
Germany

Phone: +49 61 72 608-0
Fax: +49 6172 608-5786

composeal@fresenius-kabi.com
www.fresenius-kabi.com