

Agilia Rolling Stands for Increased Mobility

Multi-Channel Rolling Stand

Specifically designed for the mobility of Link Agilia ensuring a safe centre of gravity.



Doubled pole to ensure well positioned centre of gravity
5 doubled castors with brakes
Robust weighted central base
Tube diameter: 34 mm
Telescopic bag hanger (170 to 220 cm) with 6 hooks for up to 6 kg in total
Can accommodate 8 Agilia pumps and Link Agilia racks up to a total of 21.8 kg

Art. N°: Z073160

Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Manufactured by Provita. Distributed by Fresenius Vial.



▶ Link Agilia range

Offers a robust column to improve the organisation of 4, 6 or 8 Agilia pumps



▶ Folding Handle

Optional accessory to facilitate handling and transport.

Art. N°: Z073117



▶ Tray

Optional accessory to support light objects, 2 kg max.

Art. N°: Z073113



▶ U-support

Optional accessory to fix light IV devices, 2 kg max.

Art. N°: Z073116



Agilia Rolling Stands for Increased Mobility

Infusion Pump Rolling Stand

Designed for the safe and stable mobility of 1, 2 or 3 infusion pumps.



- All-in-one product
- 5 doubled castors with brakes
- Robust weighted central base
- Tube diameter: 25 mm
- Telescopic bag hanger (170 to 220 cm) with 4 hooks for up to 2 kg each
- Integrated handle
- Can accommodate up to 3 Agilia pumps (or up to 6.8 kg)
- Integrated handle

Art. N°: Z073150

Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Manufactured by Provita. Distributed by Fresenius Vial.

From General Wards to Intensive Care units, the Agilia family offers a complete range of infusion pumps, disposables and accessories to cover all major IV drug administration needs.

▶ Agilia pumps

Injectomat Agilia



Volumat Agilia



▶ Agilia MC pumps

Injectomat MC Agilia



Volumat MC Agilia



▶ Agilia TIVA pumps

Injectomat TIVA Agilia



▶ Administration Sets


Volumat Lines



Agilia Rolling Stand for Increased Mobility

Twin Link Rolling Stand


Specifically designed
for the mobility of 2 Link.



Robust weighted central base to ensure stability
4 doubled antistatic castors with brakes
Tubes diameters: 38 mm
Device weight: 46 kg
Telescopic bag hanger (184 to 220 cm)
with 6 hooks
2 additional hooks on the Twin Link structure
Up to 2 Link Agilia and 16 Agilia pumps
Max devices weight on each poles: 25 kg (50 kg total)
Max container weight on each hook: 2 kg (16 kg total)

Art. N°: Z073170

Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Manufactured by Provita. Distributed by Fresenius Vial.



▶ Link Agilia range

Offers a robust column to
improve the organisation of
4, 6 or 8 Agilia pump

