

As of 12/30/24

[Product Certificates](#)

[X4R5720 Product Insert \(Bioflex - BI\)](#)

[X4R5720 Product Insert \(Sepacell - SE\)](#)

A new leukoreduction filter will be incorporated into the Alyx X4R5720 kit. Visual indicators (BI or SE) will be added to the Master Carton box and to the Red Cell Storage Container label to identify the kits produced with Bioflex RC XL (BI) filters versus Sepacell (SE) filters. Please reference the list below for finished good lots produced with Bioflex (BI) or Sepacell (SE) filters.

For further questions, please reach out to your Fresenius Kabi representative.

<u>Manufacturing Date</u>	<u>Code</u>	<u>Lot</u>	<u>Filter</u>	<u>Label Indicator</u>
4/29/2024	X4R5720	FN24D29011	Bioflex	None
4/30/2024	X4R5720	FN24D30019	Bioflex	None
8/13/2024	X4R5720	FN24H12018	Bioflex	BI
8/15/2024	X4R5720	FN24H14014	Bioflex	BI
8/20/2024	X4R5720	FN24H20037	Sepacell	SE
8/22/2024	X4R5720	FN24H22017	Sepacell	SE
8/23/2024	X4R5720	FN24H23015	Sepacell	SE
8/26/2024	X4R5720	FN24H26018	Sepacell	SE
8/27/2024	X4R5720	FN24H27016	Sepacell	SE
8/28/2024	X4R5720	FN24H28014	Sepacell	SE
8/29/2024	X4R5720	FN24H29012	Sepacell	SE
8/30/2024	X4R5720	FN24H30010	Sepacell	SE
8/31/2024	X4R5720	FN24H31018	Sepacell	SE
9/3/2024	X4R5720	FN24I03015	Sepacell	SE
9/4/2024	X4R5720	FN24I04013	Sepacell	SE
9/5/2024	X4R5720	FN24I05010	Sepacell	SE
9/6/2024	X4R5720	FN24I06018	Sepacell	SE
9/7/2024	X4R5720	FN24I07016	Sepacell	SE
9/10/2024	X4R5720	FN24I09020	Sepacell	SE
9/11/2024	X4R5720	FN24I11018	Sepacell	SE
9/13/2024	X4R5720	FN24I13014	Sepacell	SE
9/14/2024	X4R5720	FN24I14012	Sepacell	SE
9/16/2024	X4R5720	FN24I16017	Sepacell	SE
9/17/2024	X4R5720	FN24I17031	Sepacell	SE
9/18/2024	X4R5720	FN24I18021	Sepacell	SE
9/19/2024	X4R5720	FN24I19011	Sepacell	SE
9/20/2024	X4R5720	FN24I20019	Sepacell	SE
9/23/2024	X4R5720	FN24I23013	Sepacell	SE
9/24/2024	X4R5720	FN24I24011	Sepacell	SE
9/25/2024	X4R5720	FN24I25018	Sepacell	SE
9/26/2024	X4R5720	FN24I26016	Sepacell	SE
9/28/2024	X4R5720	FN24I27014	Sepacell	SE
9/29/2024	X4R5720	FN24I29010	Sepacell	SE
10/1/2024	X4R5720	FN24J01016	Sepacell	SE
10/2/2024	X4R5720	FN24J02014	Sepacell	SE
10/3/2024	X4R5720	FN24J03012	Sepacell	SE
10/4/2024	X4R5720	FN24J04010	Sepacell	SE
10/5/2024	X4R5720	FN24J05017	Sepacell	SE
10/7/2024	X4R5720	FN24J07013	Sepacell	SE
10/8/2024	X4R5720	FN24J08011	Sepacell	SE
10/9/2024	X4R5720	FN24J09019	Sepacell	SE
10/17/2024	X4R5720	FN24J16014	Sepacell	SE

10/19/2024	X4R5720	FN24J17012	Sepacell	SE
10/24/2024	X4R5720	FN24J24018	Sepacell	SE
10/25/2024	X4R5720	FN24J25015	Sepacell	SE
10/28/2024	X4R5720	FN24J28019	Sepacell	SE
10/29/2024	X4R5720	FN24J29017	Sepacell	SE
10/30/2024	X4R5720	FN24J30015	Sepacell	SE
10/31/2024	X4R5720	FN24J31021	Sepacell	SE
11/1/2024	X4R5720	FN24K01048	Sepacell	SE
11/3/2024	X4R5720	FN24K02012	Sepacell	SE
11/4/2024	X4R5720	FN24K04018	Sepacell	SE
11/6/2024	X4R5720	FN24K06013	Sepacell	SE
11/7/2024	X4R5720	FN24K07011	Sepacell	SE
11/9/2024	X4R5720	FN24K09017	Sepacell	SE
11/11/2024	X4R5720	FN24K11013	Sepacell	SE
11/13/2024	X4R5720	FN24K13019	Bioflex	BI
11/14/2024	X4R5720	FN24K14017	Bioflex	BI
11/15/2024	X4R5720	FN24K15014	Bioflex	BI
11/18/2024	X4R5720	FN24K18018	Sepacell	SE
11/19/2024	X4R5720	FN24K19016	Sepacell	SE
11/20/2024	X4R5720	FN24K20014	Sepacell	SE
11/21/2024	X4R5720	FN24K21012	Sepacell	SE
11/24/2024	X4R5720	FN24K22010	Sepacell	SE
11/25/2024	X4R5720	FN24K25013	Sepacell	SE
11/26/2024	X4R5720	FN24K26011	Sepacell	SE
11/27/2024	X4R5720	FN24K27019	Sepacell	SE
11/30/2024	X4R5720	FN24K30013	Sepacell	SE