



**FRESENIUS
KABI**

caring for life



CompoGuard

Whole Blood Collection System

Designed to help optimize whole blood collections

The complete line of CompoGuard Whole Blood Collection systems is designed to help optimize whole blood collections in multiple collection environments.

In-line tray for optimal bag positioning and accurate weighing

Visual and audible donation prompts for an improved blood collection process

Integrated safety clamp to prevent air from entering the collection bag

Durable transport case and compact size are ideal for mobile donation settings

CompoGuard

Whole Blood Collection System

The complete line of CompoGuard Whole Blood Collection systems is designed to help optimize whole blood collections in multiple collection environments.



	CompoGuard Basic	CompoGuard Data	CompoGuard Complete
Description	<p>Basic configuration designed with key functionality and affordability in mind.</p> <ul style="list-style-type: none"> • In-line tray • Integrated safety clamp • Large display • Easy to clean soft keys 	<p>The basic model with Data configuration includes the functionality with additional offerings.</p> <ul style="list-style-type: none"> • Work station designed for operator comfort • Basic Model + DECT* option + USB port • Optional: DonationMaster Net • Data Management software system <p>*DECT=Digitally enhanced communication technology</p>	<p>The most complete offering combines accurate weighing, mixing, and sealing while also capturing data using optional DonationMaster Net Data Management software.</p> <ul style="list-style-type: none"> • Basic Model + DECT* option + Gooseneck + USB port + sealer + rechargeable battery <p>(Model shown with optional scanner)</p>
Dimensions (H x W x D)	6.5 X 8.8 X 17.5 in	6.5 X 8.8 X 17.5 in	36 X 8.8 X 17.5 in Height is with Optenna extended
Weight	7.5 lbs	10.5 lbs	18 lbs



Optional Transport Case

Dimensions (H x W x D): 11 X 9.8 X 20.8 in
Weight: 6.17 lbs

Charging function: Case allows for access to plug the AC power to the battery pack while the case is closed.

For more information please contact us at ProductSupport@fresenius-kabi.com or 1.877.779.7760.

Please refer to the instructions for use for a list of warnings and precautions associated with the device.