



**FRESENIUS
KABI**

caring for life



CompoSeal Mobilea II Portable Tube Sealer

The CompoSeal Mobilea II is a light weight, portable radio frequency (RF) sealer designed to seal medical PVC tubing in both mobile donation settings and processing environments.

The power management system uses lithium-ion battery technology to combine quality sealing performance, high battery capacity and fast charging.

High capacity lithium-ion battery

Portable, light weight design

Built for the demanding use in mobile settings

Quality seals using RF technology

CompoSeal Mobilea II

Portable Tube Sealer

The portable light weight sealer with high battery capacity

- Lithium-ion battery provides up to 1,000 seals when fully charged (validated using common sized blood bag tubing)
- Designed for portable use: the light weight power pack weighs 1.8 lbs; a robust aluminum metal housing and rubber padding improves durability in a busy environment
- Recharge an empty battery in approximately 2 hours
- Consistent quality seals
- Quick and easy operation
- Automatic adjustment of sealing time depending on tube thickness



Description	Weight	Dimensions (W x D x H)
CompSeal Mobilea II set	5.36 lbs	15.4 x 4.3 x 11.8 in
Includes:	1.8 lbs	4.3 x 2.1 x 6.9 in
Powerpack (battery and generator)	2.6 lbs (weight in use)	8.2 x 3.5 x 1.3 in
Hand-held sealer		0.7 lbs
Cable	0.1 lbs	57.4 x 0.3 x 0.3 in
Charger	0.5 lbs	2.5 x 1.5 x 4.1 in
Transport case	2.0 lbs	15.3 x 4.3 x 11.8 in

For more information please contact us at ProductSupport@fresenius-kabi.com or 1.877.779.7760.

Please refer to the instructions for use for a list of warnings and precautions associated with the device.



**FRESENIUS
KABI**

caring for life

Fresenius Kabi Canada Ltd.
165 Galaxy Boulevard, Suite 100
Toronto, ON M9W 0C8
Tel: 1.877.779.7760
www.fresenius-kabi.ca