

ENLock Adaptor

THE SAME BUT DIFFERENT



**FRESENIUS
KABI**

caring for life



IT MAY LOOK THE SAME BUT IT WILL MAKE A DIFFERENCE

THE SAME:

The ENLock adaptor functions like the old step adaptor.

- Both have the same tip dimension and the same width at the center (↔)
- Both have the same reliable compatibility with enteral funnel ports



BUT DIFFERENT:

- The ENLock adaptor prepares you for the switch to the complete ENLock system
- ENLock is the new industry-developed market standard for enteral nutrition connections, designed under the technical lead of Fresenius Kabi in collaboration with Nestlé HealthScience and Nutricia Advanced Medical Nutrition

This sticker on the boxes identifies the ENLock adaptor.

The ENLock adaptor is removable from the Luer-compatible enteral giving set during the transition period.

This prepares you for the switch to ENLock, the new enteral nutrition connection system.