

# MEDICAL DEVICES

## Agilia Holder for Ambulance

DATA SHEET



**FRESENIUS  
KABI**

caring for life

# Agilia Holder for Ambulance

## INTENDED USE

Agilia Holder for Ambulance is a medical device intended to be used in road ambulances equipped with power source (AC or DC) and a horizontal rail in order to fix one Agilia IV pump.




COMPATIBLE PUMPS	Injectomat Agilia, Injectomat MC Agilia, Injectomat TIVA Agilia syringe pumps and Volumat Agilia, Volumat MC Agilia volumetric infusion pumps.
TYPE OF AMBULANCE	Road ambulance Type B, Type C
TYPE OF SUPPORT	Horizontal rails: - Minimum width: 160 mm, - Maximum height: 45 mm - Maximum thickness: 10 mm - 10 mm distance between road ambulance wall and rail

## USE ENVIRONMENT

PRESSURE OPERATING RANGE	700 hPa to 1060 hPa 525 mmHg to 795 mmHg
STORAGE AND TRANSPORT PRESSURE	500 hPa to 1060 hPa 375 mmHg to 795 mmHg
HUMIDITY OPERATING RANGE	20% to 90% no condensation
STORAGE AND TRANSPORT HUMIDITY	10% to 90% no condensation
ALTITUDE	maximum 3000 m

## TECHNICAL SPECIFICATIONS

MAIN MATERIALS	- Aluminium for the Agilia Holder for Ambulance - Rubber for bumpers - Stainless steel for screws and axe	
H/W/D	254.5 x 265 x 187.5 mm (closed)	
WEIGHT	Approximately 1.5 Kg	
COMPLIANCE WITH CE STANDARDS	Mark according to the 93/42/CEE Medical directive EN 1789 Medical vehicles and their equipment - Road ambulances	
ACCOMPANYING TOOLS	- one hexagonal 6 pans wrench, - one 13 flat wrench, - two locking keys	
ORDER NUMBER	z073290	

Due to our policy of continuous product development as well as changes in standards, the features described are subject to change. Please contact us for current situation.

From General Wards to Intensive Care units, the Agilia family offers a complete range of infusion pumps, disposables and accessories to cover all major IV drug administration needs.

## ▶ Agilia pumps

### Injectomat Agilia



### Volumat Agilia



## ▶ Agilia MC pumps

### Injectomat MC Agilia



### Volumat MC Agilia



## ▶ Agilia TIVA pumps

### Injectomat TIVA Agilia



## ▶ Administration Sets

### Volumat Lines



Agilia is a registered trademark by Fresenius Kabi in selected countries.