

April 15, 2025

## **Compoflex Transfer Pack Container with Coupler – 600 mL** **Backorder Notification**

Dear Valued Customer,

Fresenius Kabi Canada regrets to advise that we are encountering a supply interruption on our **Compoflex Transfer Pack Container with Coupler – 600 mL**.

<b>Fresenius Kabi Product Code</b>	<b>Product Description</b>	<b>Estimated Stockout Date</b>	<b>Estimated Availability Date</b>
R6R2023	<b>Transfer Pack Container with Coupler – 600 mL</b>	April 15, 2025	75% allocation May 9, 2025

We recognize the difficulties and inconvenience that you may encounter as a result of our supply interruption, and we thank you for your patience and understanding in this matter. In the interim, please consider the use of our **Compoflex 1F Empty Container with Coupler – 1000 mL (R6R2031)** as a substitute.

If you have any questions or concerns, please do not hesitate to contact your Fresenius Kabi representative.

Sincerely,

Lucy Chen  
Sr Marketing Manager, MedTech  
Fresenius Kabi Canada